


Labeling an Enteric Outbreak Kit

Complete all information on big sticker of kit. DO NOT REMOVE. This is for public health lab use.

DO NOT USE kit if beyond expiry date.

Return expired kits to SMDHU



Enter case name:
LAST NAME, first name

Date specimen collected.

Earliest onset date of symptoms

Case health card number (if available)

Outbreak number provided by SMDHU.
2260 – 200x - ###

This sticker needs to remain here for public health lab use.
DO NOT REMOVE

These stickers will be placed on bottles inside kit by the facility.

Place the case's filled and labeled bottles in this bag. Seal shut by peeling off the blue strip. Refrigerate package (do not freeze) immediately.