

## *Instructions on how to pack an insulated hard sided cooler*



1) Approved insulated hard sided cooler



2) Frozen ice packs to cover bottom, cover with one of layer bubble wrap



3) Layer chilled cryopak blanket on bubble wrap, place thermometer probe on the blanket



4) Fold cryopak blanket over thermometer probe, close lid of cooler



5) Monitor until temperature is between 2° and 8° C



6) Once desired temperature range has been reached, place vaccine between layers of cryopak blanket, fold blanket over vaccine



7) Close lid, cooler is now ready for transportation