

Steps for instrument Cleaning and Disinfecting

1

Soak



Soak used
instruments

2

Clean



Clean using detergent,
gloves and a brush

3


Rinse



Rinse with
clean water

4

Dry



Air dry or dry with
lint free cloth

5

Immerse



Immerse in approved
disinfectant. Follow
manufacturers
instructions

6


Rinse



Rinse with
clean water

7

Dry



Air dry or dry with
lint free cloth

8

Store



Store in a dry,
clean space



**simcoe
muskoka**
DISTRICT HEALTH UNIT

Tel: 705-721-7520
Toll free: 1-877-721-7520
www.simcoemuskokahealth.org
Your Health Connection